



Calibration of Volumetric Glassware

Birudugadda Srivibhav
IIT Gandhinagar
Gandhinagar, Gujarat, India
birudugadda.srivibhav@iitgn.ac.in

Ankeshwar Ruthesha
IIT Gandhinagar
Gandhinagar, Gujarat, India
ankeshwar.ruthesha@iitgn.ac.in

Kirtankumar Patel
IIT Gandhinagar
Gandhinagar, Gujarat, India
kirtankumar.patel@iitgn.ac.in

Nikhil Kumar Lal
IIT Gandhinagar
Gandhinagar, Gujarat, India
nikhilkumar.lal@iitgn.ac.in

Abstract—This report details the calibration process for volumetric glassware, such as volumetric flasks, pipettes, and measuring cylinders. The calibration was carried out by utilizing the density of water at room temperature. Furthermore, the volume of an unidentified test tube was measured using similar calibration methods.

I. INTRODUCTION

Calibration is essential for ensuring the precision and dependability of measurements made with laboratory equipment. This is particularly critical for glassware like volumetric flasks, pipettes, and measuring cylinders, which are commonly used in titrations and other accurate chemical analyses. In this experiment, the volumetric glassware was calibrated using the known density of water at room temperature. We also determined the volume of an unidentified test tube.

II. THEORY

To calibrate volumetric glassware, the mass of the water contained in the glassware is measured and then divided by the known density of water at a given temperature. The systematic error is identified by comparing the calculated volumes with the volumes marked by the manufacturer for each piece of glassware: a graduated pipette (10 mL), a volumetric pipette (10 mL), a volumetric flask (100 mL), and a measuring cylinder (10 mL). Additionally, a test tube with an unknown volume can be evaluated for calibration purposes. The accuracy of the calibration relies on the precise reading of the meniscus and the careful handling of the equipment, which includes a weighing machine for accurate mass measurement, wash bottle for cleaning, tissues for drying, and the use of safety gear such as eye glasses and gloves. Distilled water is used to ensure consistent density measurements during the calibration process.

III. FORMULAE USED

The primary formula used for volume calculation is:

$$V = \frac{M}{\rho} \quad (1)$$

where V is the volume, M is the mass of water, and ρ is the density of water. For standard deviations, the formula used is:

$$\sigma = \sqrt{\frac{1}{N-1} \sum_{i=1}^N (x_i - \mu)^2} \quad (2)$$

where μ is the mean of the measured values, and x_i are individual readings.

The percentage error for the graduated pipette (10 mL), regular pipette (10 mL), volumetric flask (100 mL), and measuring cylinder (10 mL) can be calculated using the formula:

$$Error = \left(\frac{\text{Actual Volume} - \text{Expt Volume}}{\text{Actual Volume}} \right) \times 100\% \quad (3)$$

The percentage error for the unknown test tube can be calculated using the formula:

$$Error = \left(\frac{\mu - \text{Expt Volume}}{\mu} \right) \times 100\% \quad (4)$$

IV. APPARATUS

The following apparatus were used:

- **Volumetric Flask:** Employed for precise measurement of a fixed liquid volume, marked with a single calibration line.

- **Volumetric Pipettes:** Used to transfer accurate volumes of liquid, particularly suited for titration purposes due to their high precision.
- **Measuring Cylinder:** A graduated container for measuring liquid volumes, useful for general purposes but less precise compared to volumetric flasks.
- **Beakers:** Cylindrical vessels used for mixing and heating liquids, featuring volume graduations, though not highly accurate for measurements.
- **Analytical Balance:** A sensitive scale for measuring mass with high precision is essential for experiments requiring exact quantities.
- **Distilled Water:** Purified water free from impurities, used to ensure the accuracy of experimental results.
- **Tissue Paper:** Used to clean glassware between iterations, ensuring that no residues affect subsequent measurements.

V. PROCEDURE

- 1) All apparatus were thoroughly cleaned and dried to ensure there were no contaminants that could affect the measurements.
- 2) Each piece of dry glassware was weighed individually using an analytical balance to obtain its initial mass.
- 3) The glassware was filled with distilled water up to the calibration mark, ensuring precise measurements for each apparatus.
- 4) The weight of the glassware, now containing water, was measured again to determine the total mass.
- 5) The mass of the water was calculated by subtracting the mass of the dry glassware from the mass of the glassware filled with water.
- 6) The volume of water was then calculated using the known density of water at room temperature, which is approximately 0.997 g/ml.
- 7) This entire procedure was repeated three times for each piece of glassware to ensure accuracy and reliability in the measurements.

VI. RESULTS AND OBSERVATIONS

In this section, we present the results obtained from the experiment, including the data acquired and calculated volumes of water.

A. Data Acquired from the Experiment

The weights measured for each apparatus, both empty and with water, are shown in Table I. The consistent increase in weight with the addition of water across all instruments indicates reliable measurements.

B. Calculated Volume of Water

The calculated volumes of water for each apparatus at a temperature of 27 degrees Celsius are summarized in the below Table II.

| Apparatus | Empty (g) | Weight (g) | | |
|-------------------------|--------------|------------|---------|---------|
| | | 1 | 2 | 3 |
| Volumetric Flask | 66.295 | 165.713 | 165.610 | 165.732 |
| Measuring Cylinder | 32.768 | 42.743 | 42.649 | 42.685 |
| Graduated Pipette | 28.284 | 38.250 | 38.289 | 38.253 |
| Volumetric Pipette | 28.284 | 38.282 | 38.269 | 38.278 |
| Test Tube (with Beaker) | 40.676 | 48.867 | 48.835 | 48.826 |

TABLE I: Data Acquired from the Experiment (Weights with and without Water)

| Apparatus | Volume (ml) | | |
|-------------------------|-------------|--------|--------|
| | 1 | 2 | 3 |
| Volumetric Flask | 99.717 | 99.613 | 99.736 |
| Measuring Cylinder | 10.005 | 9.910 | 9.946 |
| Graduated Pipette | 9.995 | 10.035 | 9.998 |
| Volumetric Pipette | 10.001 | 10.015 | 10.024 |
| Test Tube (with Beaker) | 8.215 | 8.183 | 8.174 |

TABLE II: Calculated Volume of Water

C. Final Results

The final results, including the mean, standard deviation, and percentage errors for each apparatus, are summarized in Table III. These approximate measurements of the Weights indicate the overall accuracy and reliability of the measuring instruments.

| Apparatus | Error (%) | | | Mean (ml) | Standard Deviation (ml) |
|-------------------------|-----------|-------|-------|-----------|-------------------------|
| | 1 | 2 | 3 | | |
| Volumetric Flask | 0.283 | 0.387 | 0.264 | 99.688 | 0.054 |
| Measuring Cylinder | 0.05 | 0.9 | 0.56 | 9.953 | 0.039 |
| Graduated Pipette | 0.05 | 0.35 | 0.02 | 10.009 | 0.018 |
| Volumetric Pipette | 0.01 | 0.15 | 0.24 | 10.013 | 0.009 |
| Test Tube (with Beaker) | 0.30 | 0.08 | 0.19 | 8.190 | 0.017 |

TABLE III: Calculated Volume of Water

Overall, the volumetric pipette demonstrated the highest accuracy and reliability in measuring water volumes, making it the preferred choice for precise measurements in laboratory settings. The variations observed in other apparatus suggest areas for improvement in measurement techniques and instrument selection.

VII. CONCLUSION

In this experiment, we calibrated various pieces of volumetric glassware, including a volumetric flask, pipettes, and a measuring cylinder, using the density of water at room temperature. We also determined the volume of an unknown test tube by applying similar calibration techniques. The results showed that the volumetric pipette had the highest

accuracy and smallest error, making it the most reliable for precise measurements. Others also performed well, though small deviations in volume were observed.

The volumetric flask and measuring cylinder, while adequate for general use, showed slightly higher errors, indicating they are not as accurate for highly precise tasks. These variations can be due to factors like reading the meniscus inaccurately or minor inconsistencies in the apparatus.

Overall, this experiment emphasised the importance of regular calibration in maintaining accuracy in laboratory measurements. It also highlighted the differences in precision across different types of glassware, demonstrating the need to choose the appropriate tool based on the required level of accuracy for a given task.

VIII. ACKNOWLEDGMENT

We would like to express our sincere gratitude to lab staff at IIT Gandhinagar for their invaluable guidance and support throughout this experiment. Special thanks to the laboratory staff for providing the necessary equipment and ensuring a smooth experimental process. We also extend our appreciation to our group members for their teamwork and contributions. Lastly, we acknowledge IIT Gandhinagar for providing the resources and facilities to conduct this experiment.

Handwritten experimental data table for calibration of volumetric glassware. The table includes measurements for various glassware types (beaker, volumetric flask, measuring cylinder, test tube, graduated pipette) and their corresponding errors (E₁, E₂, E₃). The data is organized into two main sections, one for 'Group 5' and another for 'Group 1'. The table is signed by TA.

| Group | Item | Value | E ₁ | E ₂ | E ₃ |
|---------|-------------------------|-----------|----------------|----------------|----------------|
| Group 5 | beaker | 28.284 gm | | | |
| | Volumetric flask | 66.295 gm | 165.713 | 165.610 | 165.732 |
| | measuring cylinder | 32.762 gm | 42.743 | 42.649 | 42.685 |
| | test tube | 12.383 gm | | | |
| | test tube (with broken) | 48.867 | 48.835 | 48.826 | |
| | graduated pipette | 28.289 | 38.289 | 38.289 | 38.289 |
| | test tube | 48.867 | 48.835 | 48.826 | |
| | graduated pipette | 28.289 | 38.289 | 38.289 | 38.289 |
| | test tube | 48.867 | 48.835 | 48.826 | |
| | graduated pipette | 28.289 | 38.289 | 38.289 | 38.289 |
| Group 1 | Volumetric flask | 99.418 | 99.315 | 99.437 | |
| | measuring flask | 99.975 | 99.881 | 99.917 | |
| | test tube | 8.191 | 8.159 | 8.150 | |
| | Bulk pipette | 9.998 | 9.985 | 9.994 | |
| | graduated pipette | 9.966 | 10.005 | 9.969 | |

Fig. 1: Experimental Data (Signed by TA)